

Hiflex⁵ the perfect complement for your therapy furniture – simply modern and functional



Hiflex⁵ – wherever it is needed



The Hiflex⁵ can be moved easily in every desired position and the base frame with its small shell matches the Bionic Therapy Chairs perfectly.



The multifunctional table can be height adjusted smoothly and fast; the big tray is tiltable in two steps.



There is a practical hook made of Aluminium for clothes or bags on the inside of the column. The Hiflex⁵ can be equipped with a useful container or...



..., if desired, with a big newspaper holder. Two separate trays offer much space for medical equipment or personal belongings during treatment.

With the simple and continuous height adjustment by foot pedal the Hiflex⁵ meets your patients' wishes. Additionally, the table top is tiltable up to 30° in two steps.

Four smooth-running wheels as well as a practical handle allow an easy movement and positioning of the Hiflex⁵.

At the inner side of the column there is a stable hook made of Aluminium for storage of bags or clothes. Alternatively, the Hiflex⁵ is equipped with a container or a newspaper holder and therefore offers enough space.

The Hiflex⁵ is available in light grey (RAL 7035) to match the Bionic Therapy Chairs perfectly.

Technical Data:

Dimensions

Total width	40 cm
Total length	78 cm
Height adjustment	74 – 97 cm

Removable Trays

Big tray	55 x 39 cm
Small tray	15 x 39 cm

Safety

Safe maximum load	15 kg
-------------------	-------

Optional Accessories

Container, H x L x W	22 x 20 x 30 cm
Newspaper holder, H x L x W	22 x 13 x 31 cm

Available colours

Light grey	RAL 7035
------------	----------

