

Order code	Description	Packaging unit	Connector material
AV-EXT2	Connection aid for double lumen catheters	50 units / carton	PC, PP
AV-EXT2-Y30	Connection aid for fistula needles	80 units / carton	PC, PP
CK311	Compression clamp	1 pair	
CLAMP-WK	Vas-Clamp compression strap with hook-and-loop fastener, sterilisable, white	1 pair	
CON-FF1	Circulation connector, Luer Lock f/f, 5 cm	100 units / carton	PC
CON-MM1	Circulation connector, Luer Lock m/m, 5 cm	50 units / carton	PC
DKS12	Dilator, 2-stage, from 10 to 12 F, 200 mm	1 unit	
DKS14	Dilator, 1-stage, 14 F, 200 mm	1 unit	
DT-A-S	Demers-Katheter® Tunneler, atraumatic	1 unit	
DT-S-S	Demers-Katheter® Tunneler, sharp	1 unit	
DPX-16	Introducer set for Demers-Katheter®	1 unit	
EXT150LL	Extension line with Luer Lock, 150 cm	25 units / carton	PC
G700	Guide wire with straight tip, 0,95 mm x 700 mm	1 unit	
G700J	Guide wire with J-tip, 0,95 mm x 700 mm	1 unit	
IKLL-ISO	Injection cap LL, sterile, latex-free membrane	100 units / carton	PP
INF150LL	Infusion system with Luer Lock, 150 cm	25 units / carton	PC
KAB21	Sharps container 2.1 l	100 units / carton	
KAB75	Sharps container 7.5 l	40 units / carton	
KPA06006-4	Connection and disconnection set for single lumen catheter (customized sets on request)	20 units / carton	
KPA06004-4	Connection and disconnection set for double lumen catheter (customized sets on request)	15 units / carton	
K2565	Plastic scissor clamp	50 units / carton	
LWS5ROT	Replacement connector for Demers-catheter® including tube clamp	1 unit	PVDF
MSK-G	MultiClamp	50 units / carton	
SNA3	Single-Needle-Adapter with Twist-Lock connector and clamp, 10 cm tube	50 units / carton	PP, PVC
SN57	Seldinger puncture Needle	25 units / carton	
SY100LL	Disposable syringe 100 ml	25 units / carton	PP
UVL	Universal cap, Luer Lock, red	100 units / carton	ABS
YLL	Y-Adapter with LL, sterile, without safety cap	50 units / carton	PC

Bionic-Accessories

secure, safe and effective

- for a secure dialysis

Universal cap, circulation connector, infusion systems, Single-Needle-Adapter, scissor clamps

- for a safe catheter placement

Introducer set for Demers-Katheter®, tunneler for percutaneous insertion, puncture needles, dilators

- for an effective connection and disconnection

Sets for single and double lumen catheters, customized sets on request



Bionic Medizintechnik GmbH ... a JMS-corporation • Max-Planck-Straße 21 • D-61381 Friedrichsdorf/Germany
Phone +49 (0) 6172-7576-0 • Fax +49 (0) 6172-7576-10 • info@bionic-jms.com • www.bionic-jms.com

Right for technical modifications reserved.
B-ALLZUE/11 2014-12



Bionic-Accessories



AV-EXT2
Connection aid for dual lumen catheters, each 10 cm long. The clamps in red and blue assure safe and fast clamping.



AV-EXT2-Y30
Connection aid for fistula needles, each 30 cm long. Arterial line with y-fitting to add saline solution at the same time.



CK311
The compression clamp for an effective compression of puncture channels after dialysis.



CLAMP-WK
The Vas-Clamp compression strap for dialysis patients is the ideal aid for blood-stanching.



CON-FF1
Circulation connector with 5 cm length with Luer Lock connectors female/female.



CON-MM1
Circulation connector with 5 cm length with Luer Lock connectors male/male.



DKS12 / DKS14
One or two-stage dilators available in different dimensions.



DT-A-S / DT-S-S
Tuner with sheath for percutaneous insertion of Demers-catheters® with atraumatic or sharp tip.



DPX-16
Introducer set for percutaneous insertion of Demers-catheters®.



EXT150LL
Extension line with Luer Lock connectors, 150 cm.



G700 / G700J
Guide wire with straight tip or j-tip.



IKLL-ISO
Latex-free injection cap with self-sealing injection membrane for Luer Lock connectors. Available in yellow.



INF150LL
Infusion system for infusion bags, 150 cm.



KAB21 / KAB75
Sharps container with a volume of 2.1 or 7.5 liter.



KPA06006-4 / KPS06004-4
Connection and disconnection set for single or double lumen catheter contains different components, e.g. drape.



K2565
The plastic scissor clamp for easy clamping of tubes of up to 7 mm diameter.



LWS5ROT
Replacement connector for the original Demers-catheter® including a white clamp.



MSK-G
MultiClamp for gentle clamping and securing of tubes as well as unscrewing of Luer-Lock connections.




SNA3
Y-Adapter with two connection tubes, each 10 cm long. The colour coded clamps assure fast and safe clamping.



SN57
Puncture needle for percutaneous insertion according Seldinger.



SY100LL
Disposable syringe 100 ml with Luer Lock connector.



UVL
Universal cap for Luer Lock connectors. Available in red.



YLL
Suitable Y-Adapter with dual connection for Luer Lock connectors.