

MultiClamp

- > gentle clamping of tubings
- > easy tubing fixation
- > simple loosening of Luer-Lock connections
- > flexible carrier for tubing systems




Versatile usage - Bionic MultiClamp



The Bionic MultiClamp is used for a simple and gentle clamping of different tubings.



Tubings can also be fixed safely and easily with this multifunctional scissor clamp.



Firmly closed Luer-Lock connections are loosened fast and effortless.



The partly opened handle is used for individual mounting of the MultiClamp and as a support for various tubing systems.

Order code	Description	Packaging unit
MSK-G	Multifunctional scissor clamp MultiClamp, green	50 units / carton

