

**For your patient's comfort
choose the PremiumLine Oncology Chair**



PremiumLine

bionic

Premium comfort, functionality and privacy for your oncology patients



The neck pillow of the large Premium headrest supports a comfortable head position, while the tiltable side segments protect your patient's privacy.



An unique armrest design creates a high class atmosphere. The knee cushion, which can be moved to the desired position, offers a comfortable leg support.



The telescopic footrest can be adjusted electrically to fit very small as well as tall patients and in addition it can be folded in or out as requested.



An activation of the shock position by foot pedal plus the high speed technology save valuable time in an emergency case and support the patient's safety significantly.



The remote control provides greater flexibility as all desired chair positions can be adjusted easily by the medical team or the patient directly.



The central brake safely locks all wheels simultaneously. While unlocked the chair can be moved easily with smooth going 10 cm wheels.

Within the oncology clinic the ambience is important, helping the patients to feel more relaxed.

The new Bionic PremiumLine Therapy Chair has been designed to meet the patient's needs by creating a comfortable, supportive and relaxing environment whilst not compromising on clinical functionality or safety.

The exclusive design of the PremiumLine Oncology Chair underlines the high quality standards that are at the core of your oncology clinic.

The comfortable upholstery creates a pleasant environment and feeling of well being for your patients. Five electric motors offer simple and smooth adjustment of all chair positions further adding to patient relaxation whilst receiving treatment.

The innovative design of both the head support and the wide armrests offer increased privacy and comfort.

The large electrically adjustable footrest provides a more secure position for the patient during treatment and can be easily folded in when not in use.

The smooth electrical height adjustment is user friendly and makes the working environment of the clinical team more ergonomic. For maximum safety the high speed shock position can be activated by a foot pedal.

A variety of accessories can be mounted on the chair to give practical space saving solutions.

The optional LED lamp can be switched on/off with the remote control and provides a pleasant light which allows your patients to read or work during treatment.

A wide range of upholstery colours are available to match your clinic's interior design concept.

The PremiumLine Therapy Chair with its unique armrest design and the headrest with foldable side segments creates a comfortable environment and improves the feeling of relaxation and privacy within your oncology suite.

Accessories to meet individual requirements



The LED lamp provides an individual reading light and helps to create a more pleasant reading and working atmosphere.



The IV pole can be mounted easily and can be a space saving alternative to a free standing version.



The side table is also space saving and offers room for personal belongings or medical equipment.

Crème (Standard colour)

Pistazio

Dolphin

Safran

Kirsche

Ocean

Available colours: Up to 35 upholstery colours are available. Please contact us!

Technical Data:

Dimensions

Seat segment width	65 cm
Chair width, incl. armrests	95 cm
Total length	220 cm
Total weight	96 kg
Ground clearance	6 cm

Electrical specifications

Protection class	2
Protection category	IPX 4
In compliance with	EN 60601
Motors	5

Adjustment features

Back segment angle	75 – 0°
Seat segment angle	30 – 5°
Leg segment angle	5 – 70°
Footrest	15 cm
Height adjustment	60 – 80 cm

Safety

Safe maximum load	180 kg
Angle of shock position	– 12°
Time to shock position	10 sec



Bionic Medizintechnik GmbH ... a JMS-corporation · Max-Planck-Straße 21 · D-61381 Friedrichsdorf/Germany
Phone: +49(0)6172-75760 · Fax: +49(0)6172-757610 · info@bionic-jms.com · www.bionic-jms.com